

=> d his

(FILE 'HOME' ENTERED AT 10:01:53 ON 13 JUL 2001)  
FILE 'CA' ENTERED AT 10:02:00 ON 13 JUL 2001

L1 0 S WILHEMI M?/AU AND 1971/PY  
L2 1 S WILHELM I M?/AU AND 1971/PY  
L3 4 S WILHELM I M?/AU AND AUTOMATIC/TI  
L4 912 S AUTOMAT?(4A)MASS SPECTRO?  
L5 41 S L4 AND VACUUM  
L6 6 S L4 AND LOCK?  
L7 12 S L4 AND SAMPLE(3A) (LOAD? OR CASSETTE)  
L8 860 S L4 NOT L5-7  
L9 518 S L8 NOT(AIR OR CHROMATOG?)  
L10 8 S L9 AND NOVEL/TI  
L11 612 S AUTOMAT?(4A)MASS SPECTRO?/TI, IT, ST  
L12 1 S L10 AND L11  
L13 335 S L9 AND L11  
L14 10 S L13 AND SOLID  
L15 65 S L3, L5-7, L12, L14  
L16 44 S L15 NOT PY>1994

=> d 116 bib, ab 1-44

L16 ANSWER 7 OF 44 CA COPYRIGHT 2001 ACS  
AN 114:114144 CA  
TI A precision carbon-14 accelerator mass spectrometer  
AU Purser, K. H.; Smick, T. H.; Purser, R. K.  
CS US-AMS Corp., Topsfield, MA, 01983, USA  
SO Nucl. Instrum. Methods Phys. Res., Sect. B (1990), B52(3-4), 263-8  
AB Details will be presented of a high throughput, automated accelerator mass spectrometer. The system, which has been designed to provide  $^{14}\text{C}/^{12}\text{C}$  and  $^{13}\text{C}/^{12}\text{C}$  isotopic ratio information for milligram samples of graphite, can provide data with accuracies comparable to that which can be achieved using the techniques of precision beta-decay counting. As many as sixty bar-coded samples/stds. are loaded into the system at one time, which can operate completely unattended and collect and analyze data for at least one day. By using dual sources, the throughput of the instrument will be close to 5000 samples per yr.  $^{14}\text{C}$  backgrounds are expected to be well below that of sample contamination, and should allow measurements to be made on samples having an age of more than 60,000 yr.

L16 ANSWER 10 OF 44 CA COPYRIGHT 2001 ACS  
AN 106:27085 CA  
TI Device for automatic mass spectrometric analysis  
IN Yano, Masayoshi  
PA Hitachi, Ltd., Japan  
SO Jpn. Kokai Tokkyo Koho, 4 pp.  
PI JP 61058153 A2 19860325 JP 1984-178343 19840829  
AB The title device consists of several probes installed on a rotatable disk and perpendicular to it in a vacuum chamber. Each probe is automatically moved into an ionization chamber where its preheated sample is ionized for mass spectrometry and then removed for the next probe sample anal. Each probe anal. takes 1-2 h. There is no introduction of air into the system during the anal., and therefore there is no change in the spectrogram backgrounds caused by moisture in the air, and no counterflow of vacuum pump oil.

L16 ANSWER 20 OF 44 CA COPYRIGHT 2001 ACS

AN 91:101547 CA  
TI A fully automated mass spectrometer for the analysis of organic solids  
AU Hillig, Heinrich; Kueper, Hendrik; Riepe, Wolfgang; Ritter, Hans Peter  
CS Inst. Spektrochem. Angew. Spektrosk., Dortmund, D-4600/1, Fed. Rep. Ger.  
SO Anal. Chim. Acta (1979), 112(2), 123-32  
AB The automation of a mass spectrometer-computer system to process up to 30 samples without attention after sample loading is described. An automatic sample changer introduces the samples successively into the ion source by means of a direct inlet probe. A process control unit detcs. the operation sequence. Computer programs are available for the hardware support, system supervision, and evaluation of the spectrometer signals. The most essential precondition for automation, automatic evapn. of the sample material by electronic control of the total ion current, is satisfactory. The system operates routinely overnight in an industrial lab., so that day work can be devoted to different anal. problems. The cost of routine analyses is halved.

L16 ANSWER 25 OF 44 CA COPYRIGHT 2001 ACS

AN 89:96355 CA  
TI Experiments with an automatic mass spectrometer in the isotopic analysis of nuclear fuels  
AU Koch, L.; Brandalise, B.; Rijkeboer, C.; Romkowski, M.; Wilhelmi,\*\*\* M.; Brachmann, K.; Heinen, G.  
CS Eur. Inst. Transuranium Elem., Leopoldshafen, Ger.  
SO Adv. Mass Spectrom. (1978), 7B, 1052-61  
AB The title automatic mass spectrometer is described and its performance is discussed with 9 refs.

L16 ANSWER 30 OF 44 CA COPYRIGHT 2001 ACS

AN 81:165750 CA  
TI Automatic mass-spectrometric analysis. Preliminary report on development of a novel mass-spectrometric system for biomedical applications  
AU Dreyer, W. J.; Kuppermann, A.; Boettger, H. G.; Giffin, C. E.; Norris, D. D.; Grotch, S. L.; Theard, L. P.  
CS Jet Propul. Lab., California Inst. Technol., Pasadena, Calif., USA  
SO Clin. Chem. (1974), 20(8), 998-1002  
AB A mass spectrograph, coupled to automatic sample prepn. devices and ion species simultaneously, was described. The app. could permit simultaneous multicomponent anal. of appropriately prepd. samples and could detect as little as 10-15 g of a single component. This approach offered significant advantages over other methods, including conventional mass spectrometry.

L16 ANSWER 32 OF 44 CA COPYRIGHT 2001 ACS

AN 78:92116 CA  
TI Automatic analysis of uranium and plutonium in solutions  
AU Von Baeckmann, A.; Neuber, J.; Wilhelmi, M.; Koch, L.  
CS Nucl. Res. Cent., Karlsruhe, Ger.  
SO Anal. Methods Nuclear Fuel Cycle, Proc. Symp. (1972), Meeting Date 1971, 329-41 Publisher: IAEA, Vienna, Austria.  
AB A system based on x-ray fluorescence anal. and mass spectrometry for detn. of U and Pu in solns., e.g. in reprocessing plants, consists of automatic sample prepn. stage in which exactly weighed amts. of the sample solns. are mixed with known amts. of Th, and aliquots of the samples are taken both for mass spectrometric isotopic diln. anal. and for detn. of the U and Pu concns. by x-ray fluorescence spectroscopy; automated x-ray fluorescence spectrometer for the rapid detn. of the U and Pu concns. in the prepd. sample by comparison with the Th internal std.; central lab. with an automated mass spectrometer for subsequent measurement of the ratio between

U and Pu isotopes and Th to prevent the samples from undergoing changes during storage, they are dried in small Al capsules in the sample prepn. stage. The design of the system is described and its advantages, esp. for purposes of nuclear safeguards, are outlined.

=> log y

STN INTERNATIONAL LOGOFF AT 10:20:53 ON 13 JUL 2001